

Outdoor



Stannah is a family-run British company which has been making pioneering products since 1860. We are the world's leading stairlift manufacturer with an unrivalled reputation for quality, integrity and service. We uphold this reputation by applying recognised standards to every product we make.

Outdoor

The Outdoor stairlift for straight stairs is specially designed for the great outdoors. Whether you have a long flight of outside stairs or just a few steps going down to the basement or up to the garden, installing the Outdoor stairlift is a comfortable, safe and reliable way to ensure external steps are never a problem again.

The chair

The Outdoor stairlift is specifically designed to work in a variety of weather conditions, be it rain, sunshine or snow. As the weather is increasingly unpredictable and can quickly change, the Outdoor stairlift has been tested to ensure it operates smoothly from -10°C up to 40°C.

The Outdoor stairlift is operated by a light-touch joystick that can be positioned on either arm and has a manual swivel seat to aid you getting off the chair at the top of the stairs. It also has safety sensors that automatically stop the stairlift if they detect an obstruction on the stairs, and a seat, footrest and arms that fold away to leave the stairs clear for other users.

It is also powered by a battery that continually charges automatically when the stairlift is parked, which means it's always ready for use, even in the event of a power cut.

We offer a protective cover to keep your stairlift clean when not in use. We can also install additional call stations so you can bring the lift to you or send it up the stairs to someone else.

The rail

The Outdoor stairlift runs along a strong rail with a slim line design and it is fitted to the stairs and not to the wall. Due to the conditions in which it will be used, the Outdoor stairlift can only be installed into solid concrete, solid wood or metal staircases.

The service

We pride ourselves on the reliability of our stairlifts. To give you absolute peace of mind, they all come with a full guarantee, a free first service, and support from our nationwide teams of locally-based engineers, who are on standby 24-hours a day, 365 days a year to respond quickly in the unlikely event that you need assistance.



1. To use the Outdoor stairlift, unfold the footrest, the seat and the arms.



2. Once unfolded, make yourself comfortable on the seat. Use the safety belt for extra security if you choose.



3. Gentle pressure on the directional controls positioned on the arm is all it takes to make your stairlift move smoothly up or down the stairs.



4. When you reach the top of the stairs, you can swivel the chair round to face the landing. To ensure your safety, the chair will only ever swivel to face up the stairs.



5. After dismounting, you can store the chair away by simply lifting up the arms, seat and footrest. Or, you can leave the chair in the swivel position so it acts as a safety barrier to prevent anyone falling down the stairs.



6. When not in use, the stairlift can be folded flat against the wall so the staircase can be used by others. A protective cover is also provided.



Directional controls

can be positioned on either arm and is a light-touch joystick that is easy to use.



Manual swivel seat

is operated with easy-to-use levers on either side of the chair, allowing you to get off the chair with maximum ease.



Folding seat

can be folded away when not in use to keep the stairs clear for other users.

Self-charging battery

powers the stairlift and ensures a smooth ride with a gentle start and stop. The Outdoor stairlift charges automatically when it is parked, so it is ready for use even in the event of a power cut.

Safety sensors

on the chair stop the stairlift automatically if it meets an obstruction on the stairs.



Safe and easy alighting



Optional call station



Protective cover



Padded arms

can be folded away to leave the stairlift flat against the wall.

On/Off switch

Keyswitch

can be used to completely immobilise the stairlift, ideal if children are around.

Seatbelt

to keep you secure in the chair.



Slim strong rail

is fitted to the stairs and not to the wall. The rail can be fitted to solid concrete, solid wood or metal staircases.

Footrest

folds away to clear the stairs for other users.



Stannah contact details:
Watt Close
East Portway
Andover
Hampshire
SP10 3SD

Tel: 0800 715358
Fax: 01264 338043

www.stannahstairlifts.co.uk

Our company policy is one of continuing product improvement. We therefore reserve the right to alter the product specifications and to change and improve design at any time. Full operating instructions are supplied for all products. Stannah is a registered trademark. The CE mark indicates that this product complies with the applicable European directives which relate to health, safety, environment and consumer protection. For design, manufacture and installation, Stannah Stairlifts Limited were the first stairlift manufacturer to obtain full qualification to British Standard BS5776 (1996). An example of this product has been independently tested and approved.

