




## **SPECIFICATION**


The Poolpod Platform Lift provides dignified and independent access for swimmers who find it hard to get in and out of swimming pools.

- Independently operated by way of electronic wristband
- User can choose to either remain standing or sit using the custom submersible wheelchair
- Maximum working load 250Kg
- Powder coated 316 grade stainless steel framework
- Powered by a 24V rechargeable battery
- Semi-portable transport trolleys available to allow the Poolpod to be moved away from the poolside
- Takes just 20 Seconds to lower into the pool
- Lowers to a depth of 0.9m
- Installation requires just two anchors at the rear of the Pod
- Suitable for retrofit or new build pools
- Compatable with most pool types
  - A. Total Poolpod Length
  - B. Total Poolpod Height
  - C. Total Poolpod Width
  - D. Maximum Platform Depth Into pool












**Poolpod Platform Lift - PP1**Specification





## **SPECIFICATION**

Submersible wheelchair for use in swimming pools and wet leisure activities.

- Framework 316 grade stainless steel
- Hand operated pivot brake system
- 20" plastic rear wheels and push rails
- Pneumatic tyres
- Option of sealed stainless steel or plastic wheel bearings
- Quick dry fabric seating
- Chest strap for added lumber support.







- A. Total Wheelchair Length
- 3. Total Wheelchair Width
- C. Seating Width
- D. Foot Rest Width
- E. Distance Between Front Castors
- F. Wheel Base
- G. Seat Height
- H. Back Rest Height
- I. Seat Angle



**Submersible Wheelchair** Specification