

Junior Folding Lift

Handles all conventional lifting situations then compacts for storage under the bed or behind the door. The Junior is the ideal combination of function and convenience.







Use the ergonomically designed hand control.




Or the secondary set of controls with emergency stop mounted on the control box. Control box also has a built-in cycle counter.



Charge the battery on the lift or with an optional wall mounted charger.

AVAILABLE SLINGS



HAMMOCK SLING

For most transfers. Head-rest incorporated. Suitable for most individuals. Excellent for bathing.



QUICK-FIT SLING

For general transfers.
Perfectly suited for stocky
or obese patients with
strong hips or large thighs.



HYGIENIC SLING

Ideal for transferring to the toilet. It also permits an easy way to change, dress or undress the individual.

	SMALL	MEDIUM	LARGE
Suggested weight	45-100 lbs. (20-46 kg)	100-200 lbs. (45-90 kg)	200-600 lbs. (90-272 kg)
Suggested height	4'-5'11" (120-180 cm)	5'-5'11" (152-180 cm)	6' + (181 cm +)
Hammock	THA-S	THA-M	THA-L
Hammock 6	THA6-S	THA6-M	THA6-L
Quick-fit	TIR-S	TIR-M	TIR-L
Hygienic	THY-S	THY-M	THY-L



BHM Medical Inc. 2001 Tanguay Street Magog, Quebec, Canada J1X 5Y5 Toll free: 1-800-868-0441 Tel.: (819) 868-0441 Fax: (819) 868-2249 Internet: www.bhm-medical.com

Internet: www.bhm-medical.com E-mail: bhm@bhm-medical.com



