



HOME / PRODUCTS / MOLIFT-PASSIVE-LIFTERS / MOLIFT SMART 150 TRAVEL CASE

Molift Smart 150 Travel Case

Ensure safe transport for your Molift Smart 150 patient lifter.

Molift Products

Molift Stand-up Lifts

Molift Passive Patient Lifts

Molift Smart 150
Molift Smart 150 Travel Case
Molift Mover 205
Molift Partner 255
Molift Mover 300
Accessories

Molift Ceiling Hoist Systems

Molift Slings

Molift Stretchers



Easy Travelling

Travelling with the Molift Smart 150 is easy and convenient with this uniquely designed accessory. Produced from robust reinforced plastics, the Molift Smart 150 Travel Case protects your patient lifter from damage during travel or transit.

The Molift Smart 150 Travel Case safely stows your folded Molift Smart 150, battery, charger, and slings to ensure protection on any trip. Whether travelling by plane, train, or car, the Smart 150 Travel Case makes sure your patient lifter arrives at your destination safely and securely.

Equipped with lightweight durable rear wheels for easy carriage, the Molift Smart 150 Travel Case is easily transported to and from any location.

Portable protection for your peace of mind.

The Molift Smart 150 Travel Case is the ideal accessory for any caregiver, homecare nurse, or patient who enjoys frequent travel to family or friends. The Molift Smart 150 Travel Case provides protection for your Smart 150, which means your patient lifter can accompany you everywhere and to every destination.

Features

- Compact, lightweight and portable
- Easy-to-use wheel and move
- Provides optimum protection
- Easily stows your folded Molift Smart 150
- Supplied with lock and key for added security, TSA Approved for plane flight



Specifications

- Color: Black
- Product weight: 30lbs (This item is not 11 lbs)
- Materials: Durable reinforced plastics

Fitting the Molift Smart in the hard case

- Fold the lifter completely as described in the lifter's operator manual.
- Remember to put the suspension in the hook and the cord around the legs and castors.
- Pick the lifter up and place it into the hard case upside down aligning the wheels of the lifter base and the hard case wheels.